



Xenon light source



Microfluidic Reagent Rotor

7 minutes


Sample-to-answer in
3 simple steps and
approximately 7 minutes

Chemistry Analyzer
Celercare® cV5

MNCHIP - Clinical Chemistry +



Celercare cV5 Reagent Discs



	Health Checking Profile	Health Checking Plus Profile	Comprehensive Profile	Triple tests Profile (3)**	Preanesthetic Panel	Preanesthetic Panel Plus	Liver & Kidney Profile	Liver Profile	Kidney Profile	Electrolytes	Critical Care Panel	Avian & Reptile Panel	Large Animal Diagnostics	Diabetes Panel	Pancreatitis Profile	Equine Profile	Blood Gas Profile	Ammonia Test Profile	TBA Profile
TP	●	●	●	●	●	●	●	●			●	●		●	●				
ALB	●	●	●	●	●		●	●	●		●	●		●	●				
GLO*	●	●	●		●		●	●			●	●		●	●				
ALB/GLO*	●	●	●		●		●	●			●	●		●	●				
ALT	●	●	●	●	●	●	●	●	●										
AST			●		●	●	●	●			●	●						●	
AST/ALT*			●		●	●	●	●											
GGT			●			●	●	●					●					●	
ALP	●	●	●	●	●	●		●				●						●	
TBIL	●	●	●	●			●	●										●	
DBIL								●											
IBIL*								●											
TBA		●	●					●			●								●
CRE	●	●	●	●	●	●	●		●	●								●	●
BUN	●	●	●	●	●	●	●		●	●		●						●	●
BUN/CRE*	●	●	●		●	●	●		●	●								●	●
UA											●								
CK	●		●	●		●					●	●						●	
LDH						●													
AMY	●	●	●	●										●					
LPS														●					
GLU	●	●	●	●	●	●	●			●	●		●	●	●	●			
FRU													●						
TG			●										●						
CHOL	●	●	●	●									●						
K ⁺		●							●	●	●	●						●	●
Na ⁺		●								●	●	●						●	●
Na ⁺ /K ⁺		●									●	●	●					●	●
Cl ⁻										●	●	●						●	●
Ca	●	●	●	●					●	●	●	●		●	●	●			
P	●	●	●	●					●	●	●	●		●	●				
Ca x P*			●																
Mg			●						●			●							
tCO ₂			●						●	●	●						●	●	
NH ₃																		●	
LAC																		●	
pH																		●	

Coagulation Profile (4)*
PT APTT TT FIB

* Calculated test value
**Select three analytes from 13 test analytes for reporting